

Prevent the spread of germs and illness with the swipe of a wipe.

United 370

UNITED MD

Bleach Wipe 1:50

UNITED MD Bleach Wipe 1:50 is registered as a bactericidal effective in the presence of a 5% organic soil load against:

- *Staphylococcus aureus* Methicillin Resistant [MRSA]
- *Enterococcus faecalis* Vancomycin Resistant [VRE]
- *Streptococcus pneumoniae*
- *Streptococcus pyogenes*
- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*



Each wipe contains a 1:50 dilution of 5.25% Sodium Hypochlorite Bleach. This provides 1,000 ppm available chlorine, falling in line with CDC and OSHA recommendations.



Use on hard, non-porous surfaces including:

- Stainless Steel
- Metal
- Chrome
- Laminate
- Glass
- Sealed or Painted Wood
- Vinyl
- Linoleum
- Plastics
- Glazed Porcelain
- Glazed Ceramic

Areas of use:

- Medical, Dental and Laboratory Counters
- Exam Tables
- Carts
- Dialysis Machines
- Point of Care Equipment
- Telephones
- Sinks
- Counter Tops
- Floors
- Bathrooms
- Tubs
- Handrails
- Door Knobs
- Lockers
- Desks



Each premoistened wipe measures 8" x 10" and is packed in a fully sealed foil wrap for convenience and ease of use. Wipes are ready-to-use – no diluting or handling of concentrated bleach.

UNITED
LABORATORIES